V-Shape III

Circumference Reduction


Cellulite Reduction


Body Contouring

Skin Tightening

Breasts enlargement&Buttocks lifting








Principle



U-is a super device for shaping and removing wrinkles that integrates 5 technologies such as RF (radio frequency), VAC (negative pressure), infrared light energy, mechanical rollers, and skin floding. Through aerobic exercises such as suction, pinching and massage of connective and fat tissue, which cause deep heating of the fat cells, their surrounding connective tissue and the underlying dermal collagen fibers. Combined with 650nm near infrared and radio frequency thermal energy physical therapy, this type of efficient heating and vacuum stimulates the growth of new and better collagen and elastin which results in localized reduction in skin laxity, body volume, and an overall improvement in skin structure and texture.



Technical advantages



1 Vacuum Technology

It is used to adsorb the cellulite in a soft and uniform way, no damage to skin.



2 Multi-directional mechanical rollers

Wavelike massage speeds up blood circulation and lipolysis.



3 Unique 650nm-1200nm Infrared Light Energy

Optimize energy and accurately target the dermis, increase the safety of treatment




4 Bipolar Radio Frequency

Quickly enter the deep layer to reach a predetermined temperature, melt fat and increase collagen




5 Suction Technology

Assist the RF energy to penetrate the skin, ensure the skin to be heated evenly and safely.

















ture Health ement regimen

Application:

Fat reduction

Tighten the abdomen Shape the waist Slim limbs Breast enhancement & Buttocks Lifting

Health regimen

Meridian Treatment, Muscle Soreness Relieve,

Physical therapy
Texture improvement

Skin toning/Smooth skin, Face Contouring/Body contouring, Relieve Muscle Fatigue

- Wrinkle removal
- Neck、Forehead、eye area
 Cellulite reduction

Stretch mark, obesity lines

Unique 10 luxurious applicators

- According to the treatment site, items, demands, matching applicators. Easily achieve the purpose of perfect body figure, weight loss, skin tightening and skin texture improvement.
- Reduction of the treatment time and course of treatment to achieve the desired effect.











Handle1

Handle2

Handle3

Handle4/5/6

Handle7/8/9/10











reatment area:
Arm, Neck



Face,Eye



Instant: The size is significantly reduced after operation, Significantly Effective

One stage: The average size is reduced by 3-6 cm and the fine lines are reduced

Two weeks: Complete metabolism, weight loss, waist and abdomen tightening

After 30 days: The size reaches the expected effect, the skin is cleansed and bright, and the contours of the body appear.





